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Frozen section fluorescent single label staining experiment report (TSA)

1. Experimental principle

Tyramide signal amplification (TSA) technique is a kind of enzymatic detection method using horseradish peroxidase (HRP) to label target protein, which is similar to the conventional immunohistochemical DAB color development method. TSA technology also uses HRP labeled secondary antibodies, and also has corresponding "color rendering" steps to catalyze the addition of tyramine fluorescein substrate in the reaction system to produce an activated fluorescent substrate, which can covalently bind tyramine fluorescein and other residues on the antigen to make the sample covalently bind tyramine fluorescein. After that, the non-covalently bound primary - secondary -HRP complex was washed by thermal repair method, and the next primary -hrp secondary - antibody was repeated for the second round of incubation, so that multiple labeling could be achieved repeatedly.

2. Laboratory equipment and reagents

2.1 Laboratory equipment

Equipment	Manufacturers	Model
Frozen microtome	DAKEWE	6250
Oven	Tianjin Leibo Terry Equipment Co., Ltd	GFL125
Microwave oven	Midea	M1-L213B
Coverslips	Jiangsu Huida Medical Instruments Co.,	710510

	Ltd	
Glass microscope slides	Haimen Shenyong Experimental Equipment Factory	188109
Rocker	Servicebio	SYC-Z100
Vortex	Servicebio	MX-F
Micro-centrifuge	Servicebio	DS-S 100
Microscope	NIKON	ECLIPSE E100
Liquid Blocker PAP Pen	Gene tech	GT1001
Pipettor	Dragon	KE003068
3D Imaging System	3D HISTECH	Pannoramic SCAN

2.2 Laboratory reagents

Reagents	Manufacturers	Catlog
OCT embedding agent	Hangzhou Haoke Biotechnology Co., Ltd	HK2061
Methyl alcohol	Sinopharm Chemical Reagent Co., Ltd	67561
EDTA Antigen repair solution (9.0)	Hangzhou Haoke Biotechnology Co., Ltd	HKI0004
BSA	Hangzhou Haoke Biotechnology Co., Ltd	HKW2084
PBS solution	Hangzhou Haoke Biotechnology Co., Ltd	HK0002
Primary antibody	Hangzhou Haoke Biotechnology Co., Ltd	
Second antibody: HRP hypersensitive goat anti-rabbit secondary antibody	Hangzhou Haoke Biotechnology Co., Ltd	HKI0029
FlareXXXX	Hangzhou Haoke Biotechnology Co., Ltd	
DAPI	Hangzhou Haoke Biotechnology Co., Ltd	HKI0005

3. Experimental procedure

3.1 Frozen section fixation: The frozen slices were taken out of the refrigerator and rewarmed, dried, fixed in pure methanol for 10min, then dried at 37 °C or room temperature, and washed three times on a decolorizing shaker in PBS (PH7.4) for 5min each time.

3.2 Repair: The tissue sections were placed in a repair box filled with EDTA antigen repair buffer (PH9.0) for antigen repair in a microwave oven. Low fire for 10min, during this process should prevent excessive evaporation of buffer, do not dry. After natural cooling, the slide was placed in PBS (PH7.4) and washed by shaking on the decolorizing shaker for 3 times, 5min each time. (Repair solution and repair strength are determined according to the tissue)

3.3 Block endogenous peroxidase: The sections were incubated with 3% hydrogen peroxide at room temperature for 25min. The slides were placed in PBS (PH7.4) and washed three times on a decolorizing shaker for 5min each time.

3.4 Draw a circle around the tissue: Use a histochemical pen to draw a small circle 3-4 mm apart from the tissue along the outer outline of the tissue, and then add enough PBS to ensure that the subsequent addition of blocking serum, primary antibody, secondary antibody, and color development agent can completely cover the tissue without flowing along the sections.

3.5 Serum blocking: After the section is slightly dried, draw a circle around the tissue with a histochemical pen (to prevent the antibody from flowing away), and add 3%BSA to the circle and seal it at room temperature for more than 30 minutes.

3.6 Add primary antibody: Gently shake off the sealing solution, add PBS to the slices with a certain proportion of primary antibody, and the slices are placed flat in a wet box at 4° C for overnight incubation. (Add a small amount of water to the wet box to prevent the antibody from evaporating)

3.7 Add secondary antibody: The sections were placed in PBS (PH7.4) and washed by shaking on the decolorizing shaker for 3 times, 5min each time. After the slices were slightly dried, the tissue covered by the second antibody of the corresponding species of the first antibody was added to the ring and incubated at room temperature for 50min away from light.

3.8 TSA: The sections were placed in PBS (PH7.4) and washed by shaking on the decolorizing shaker for 3 times, 5min each time. Add the corresponding color Flare signal amplification reagent, 3-5min. Place the slide in PBS (PH7.4), shake and wash it on the decolorizing shaker for 3 times, 5min each time.

3.9 DAPI restaining nuclei: The sections were placed in PBS (PH7.4) and washed by shaking on the decolorizing table for 3 times, 5 min each time. Blot the excess PBS in the ring and then add DAPI dye solution, and incubate at room temperature for 8min away from light.

3.10 Sealing: The sections were placed in PBS (PH7.4) and washed by shaking on the decolorizing table for 3 times, 5min each time. Remove the slices and seal them with anti-fluorescence quenching tablets.

3.11 Microscopy photography: Sections were observed under a fluorescence microscope and images were collected.

4. Interpretation of results

The nuclei stained by DAPI are blue under ultraviolet excitation, and the positive expression is the color labeled by the corresponding fluorescent reagent (refer to the table below for specific fluorescence channels).

Flare Dye	Excitation wavelength	Emission wavelength
DAPI	350	420
Flare480	450	480
Flare520	490	520
Flare570	550	570
Flare620	590	620
Flare690	630	690
Flare780	750	780

5. Precautions

5.1 Don't let the slices dry out during the experiment;

5.2 Be careful of the injured tissue on the pipette suction head during experimental operation;

5.3 Pay attention to prevent peeling of frozen slices.