Western Blot 实验报告

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Western Blot Report

1. Experimental principle

Western blot(WB), is a composite immunoassay technique. The method uses SDS-PAGE technology to separate protein molecules on a gel according to their molecular weight in a biological sample, and then transfer these proteins to a solid phase membrane by means of electrical transfer. The protein on the solid phase membrane acts as an antigen and reacts with the corresponding antibody, which then reacts with an enzyme-labeled second antibody. Finally, the proteins expressed by specific target genes isolated by electrophoresis are detected by means of substrate color development or fluorescence imaging.

2. Laboratory equipment and reagents

2.1 Laboratory equipment

Equipment	Manufacturers	Model
Electronic balance	Wuxi Ying Heng Electronics Co., Ltd	YHB3003
Full wavelength enzyme labeling instrument	MolecuLar Devices	Spectramax M2
High speed refrigerated centrifuge	Haier	LX-165T2R
Laboratory pure water system	Zhiang Instrument (Shanghai) Co., Ltd	Clever-S
Magnetic stirrer	Servicebio	MS-150
Decolorizing shaker	Servicebio	DS-S100

Vertical Electrophoresis apparatus	Beijing Dongfang Rayleigh Technology Co., LTD	DYY-600C
Chemiluminescence apparatus	Clinx	ChemiScope 6100
Ultra-low temperature refrigerator	Haier	DW-86L338J
semi-dry transfer unit	ACE Biotechnology	S-TRANS FW606
Electrophoresis tank	Beijing Dongfang Rayleigh Technology Co., LTD	DYC-ZY2
Tissue Grinder	Servicebio	KZ-II
Ultrasonic cell crusher	Scientz	JY92-IIN

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网址：www.haokebio.com

邮箱：2235184086@qq.com

Centrifuge the cells, add 200 μ L RIPA per 10^6 cells (add a protease inhibitor within a few minutes before use), and oscillate. If the protein concentration needs to be increased, the volume of the total protein extraction reagent can be appropriately reduced.

3.1.2 For adherent cells:

- 1) Flush the cells with PBS 2-3 times. Drain the residue thoroughly one last time.
- 2) Add the appropriate volume of RIPA (add the protease inhibitor within a few minutes before use) in the culture plate and bottle for 3-5 minutes. During the period, the culture plate and bottle were repeatedly shaken to make the reagent fully contact with the cells.
- 3) Scrape off the cells and reagents with a cell scraper and collect them into a 1.5mL centrifuge tube.
- 4) Ice bath for 30min, during which the pipette is blown repeatedly to ensure complete cell lysis.
- 5) Centrifuge 12000g for 5min and collect the supernatant, which is the total protein solution.

3.2 Preparation of lysate from tissues :

- 1) The tissue blocks were washed with cold PBS for 2-3 times to remove the blood stain, dried with absorbent paper, cut into small pieces and placed in the grinding tube, adding a 5mm steel ball, placing the EP tube in liquid nitrogen for 2min, taking out the EP and placing it in the homogenizer for grinding (65HZ, 10s operation, 15s stop, -40°C).
- 2) When the EP tube is organized into powder, add 10 times the tissue volume of this reagent (add the protease inhibitor within a few minutes before use) and grind it in a homogenizer (65HZ, work for 10s, stop for 15s, run 20 times, -40°C).

- 3) Transfer the homogenate into 1.5mL centrifuge tube for ice bath ultrasound (200w, 15% power, ultrasound 3s, pause for 5s, run for 3min).
- 4) Centrifuge 12000g for 10min and collect the supernatant, which is the total protein solution.

3.3 Determination of protein concentration: Determination of protein concentration by BCA protein quantitative detection kit.

3.4 Denature the protein: Add 5* loading buffer to the protein solution at a ratio of 4:1, denature it in a Dry Bath Incubator at 95°C for 10 minutes, then stored at -20°C or -80°C refrigerator for later use.

3.5 SDS-PAGE Electrophoresis

- 1) Align the glass plate and put it into the clamp. Align the two glasses during operation to avoid leakage of glue. The separation glue was prepared according to the experimental arrangement. After adding TEMED, shake well immediately to fill the glue. After about 45 minutes, the water above the glue can be poured and the remaining water can be sucked up with absorbent paper.

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- 2) Mix 4% concentrated glue according to the previous method, add TEMED and shake well immediately to fill glue. Fill the remaining space with glue concentrate and insert the comb into the glue concentrate.
- 3) Add enough electrophoresis solution before loading the sample. The sample is added to the electrophoresis hole, electrophoresis. The concentrated glue voltage is 80V, and the separation glue is 120V. When the bromophenol blue has just run out, the electrophoresis can be terminated and the membrane can be transferred.

3.6 Transferring the protein from the gel to the membrane

- 1) During electrophoresis, the transferred PVDF membrane was activated in methanol for 5min, and then balanced in a incubator box containing S-TRANS anode buffer for 5min.
- 2) Place the gel that completes the electrophoresis in a tray filled with deionized water and soak for 2min.
- 3) Take two S-TRANS high-efficiency transfer gaskets and put them into the incubator box containing S-TRANS anode buffer to fully infiltrate. Put the "sandwich structure" (the gasket with anode buffer at the bottom, transfer film on top, glue on top, and cathode buffer on top) into the transfer box, and close the lid. Slide the transfer box into the transfer slot after the cancelling is fully inserted into the locking slot, 1.5A, 10min.

4) Transfer completed, remove the transfer film.

3.7 Antibody staining

- 1) The transformed film was closed for 1h with 5% skim milk (prepared with 0.5%TBST) on a decolorized shaker at room temperature.
- 2) Dilute primary antibody (TBST dissolved, phosphorylated protein using 5%BSA dissolved by TBST) and incubate overnight at 4°C.
- 3) Wash with TBST on the decolorizing shaker at room temperature five times, 5min each time.
- 4) The second antibody was diluted 10,000 times with TBST, incubated at room temperature for 120min, and washed five times with TBST on the decolorizing shaker at room temperature for 5min each time.

3.8 Chemiluminescence

Mix liquid ECL A and B in a 1:1 ratio and set aside. Take out the eluted PVDF membrane and put it on absorbent paper, slightly absorb the liquid above the membrane, put the membrane on the chemiluminescence meter shelf, add the mixed ECL luminescent liquid, let the liquid completely soak the membrane, after the reaction for 1min, use absorbent paper to absorb the excess liquid above, put it into the chemiluminescence meter, start chemiluminescence according to the preset procedure.

3.9 Analysis

Use PhotoShop software to sort out the results, and analyze the gray value of the target strip by Image J software.

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4. Results (detailed in Excel gray analysis table)

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Biotechnology Co., Ltd.